

Brabender Moisture Tester MT-CA



The Brabender Moisture Tester MT-CA uses the drying chamber method (compliant to ISO 712, ICC-Standard Nr. 110/1 etc.) and determines the loss in weight of the sample material, which results from drying.

Its new web-based MetaBridge software can be accessed by a simple browser log-in from any end device, e.g. from a PC, laptop, tablet or smartphone. From anywhere in the world, several authorized users can retrieve measurements simultaneously. Furthermore, the software allows the MT-CA to be interconnected with other USB compatible Brabender instruments.